


Anti-Stam Rabbit Polyclonal Antibody For Wb And Elisa - Q92783

Item Number: CM0030042



Introduction

Rabbit polyclonal antibody against human STAM (signal transducing adapter molecule 1). Validated for Western blot (WB) and ELISA. Suitable for human, mouse, rat samples. Target involved in cell signaling and endosomal sorting.

[Learn More](#)

Parameter	Value
Product Name	STAM Rabbit pAb
Remarks/Alias	STAM1; STAM-1; STAM
Species	Human
GeneID (Human)	8027
GeneID	8027
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 120-210 of human STAM (NP_003464.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92783
Protein Weight	59kDa
Shipping	Ice bag