

Anti-Stam2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 075886

Item Number: CM0028177



Introduction

STAM2 Rabbit Polyclonal Antibody validated for Western blot, IF/ICC, and ELISA in human, mouse, and rat. Targets signal transducing adapter molecule 2 involved in cytokine signaling and receptor down-regulation.

[Learn More](#)

Parameter	Value
Product Name	STAM2 Rabbit pAb
Remarks/Alias	Hbp; STAM2A; STAM2B; STAM2
Species	Human
Gene ID	10254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 376-525 of human STAM2 (NP_005834.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75886
Protein Weight	58kDa
Shipping	Ice bag