


Anti-Stambp Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O95630

Item Number: CM0028428



Introduction

Rabbit polyclonal antibody targeting human STAMBP (AMSH), validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Suitable for studying deubiquitination, endosomal sorting, and BMP signaling.

[Learn More](#)

Parameter	Value
Product Name	STAMBP Rabbit pAb
Remarks/Alias	AMSH; MICCAP; STAMBP
Species	Human
GeneID (human)	10617
GeneID	10617
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-270 of human STAMBP (NP_998787.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95630
Protein Weight	48kDa
Shipping	Ice bag