

Anti-Stambpl1 Polyclonal Antibody For Wb And Elisa - Q96Fj0

Item Number: CM0019612



Introduction

Rabbit polyclonal antibody against STAMBPL1 (AMSH-LP), validated for Western blot and ELISA. Detects human, mouse, and rat STAMBPL1. Ideal for ubiquitin signaling and TORC1 pathway studies.

[Learn More](#)

Parameter	Value
Product Name	STAMBPL1 Rabbit pAb
Remarks/Alias	AMSH-FP; AMSH-LP; ALMalpha; bA399O19.2; STAMBPL1
Species	Human
Gene ID (Human)	57559
Gene ID	57559
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-250 of human STAMBPL1 (NP_065850.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96FJ0
Protein Weight	50kDa
Shipping	Ice bag