

Anti-Stap2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ugk3

Item Number: CM0019654



Introduction

Rabbit polyclonal antibody targeting human STAP2 (BKS), suitable for Western blot and ELISA. Detects both human and mouse STAP2. Ideal for studying PTK6 signaling, STAT3 modulation, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	STAP2 Rabbit pAb
Remarks/Alias	BKS; STAP2
Species	Human
GeneID (Human)	55620
GeneID	55620
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-220 of human STAP2 (NP_001013863)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UGK3
Protein Weight	45kDa
Shipping	Ice bag