


Anti-Star (Stard1) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P49675

Item Number: CM0007295



Introduction

Rabbit monoclonal antibody against human StAR/STARD1, validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Targets key cholesterol transport protein in steroid hormone synthesis. Ideal for steroidogenesis and adrenal/gonadal research.

[Learn More](#)

Parameter	Value
Product Name	StAR Rabbit mAb
Remarks/Alias	STARD1; StAR
Species	Human
GeneID (Human)	6770
GeneID	6770
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human StAR (NP_000340.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49675
Protein Weight	32kDa
Shipping	Ice bag