


Anti-Star (Stard1) Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P49675

Item Number: CM0022362



Introduction

Rabbit polyclonal antibody targeting human StAR (STARD1), a mitochondrial cholesterol transporter critical for steroidogenesis. Validated for WB, IF, and ELISA; cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	StAR Rabbit pAb
Remarks/Alias	STARD1; StAR
Species	Human
Gene ID (Human)	6770
Gene ID	6770
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 64-285 of human StAR (NP_000340.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49675
Protein Weight	32kDa
Shipping	Ice bag