


Anti-Stard3 (Mln64, Cab1) Rabbit Polyclonal Antibody For Wb, Elisa - Q14849

Item Number: CM0025674



Introduction

Rabbit polyclonal antibody targeting human STARD3 (MLN64, CAB1). Validated for Western Blot and ELISA, cross-reacts with mouse and rat. Suitable for studying cholesterol transport and endosome biology.

[Learn More](#)

Parameter	Value
Product Name	STARD3 Rabbit pAb
Remarks/Alias	CAB1; es64; MLN64; STARD3
Species	Human
Gene ID (Human)	10948
Gene ID	10948
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-445 of human STARD3 (NP_006795.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14849
Protein Weight	51kDa
Shipping	Ice bag