

Anti-Stard3NI Polyclonal Antibody For Wb, Ihc-P, Elisa - 095772

Item Number: CM0020981

Introduction

STARD3NL Rabbit Polyclonal Antibody validated for WB, IHC-P, and ELISA. Detects human STARD3NL (MENTHO) protein, cross-reactive with mouse and rat. STARD3NL is a late endosome tethering protein creating ER contact sites, ideal for membrane trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	STARD3NL Rabbit pAb
Remarks/Alias	MENTHO; STARD3NL
Species	Human
Gene ID (Human)	83930
Gene ID	83930
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 185-234 of human STARD3NL (NP_114405.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95772
Protein Weight	27kDa
Shipping	Ice bag