

Anti-Stat1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P42224

Item Number: CM0002043



Introduction

STAT1 Rabbit Monoclonal Antibody validated for Western blot, immunoprecipitation, and ELISA. Reacts with human, mouse, rat STAT1. Ideal for studying interferon signaling and antiviral responses.

[Learn More](#)

Parameter	Value
Product Name	STAT1 Rabbit mAb
Remarks/Alias	CANDF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91; STAT1
Species	Human
Gene ID	6772 (Human)
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human STAT1 (P42224)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42224
Protein Weight	87kDa
Shipping	Ice bag