

Anti-Stat1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P42224

Item Number: CM0015787



Introduction

Rabbit polyclonal antibody against human STAT1 (P42224) validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying interferon signaling and STAT1 transcription factor activity.

[Learn More](#)

Parameter	Value
Product Name	STAT1 Rabbit pAb
Remarks/Alias	CANF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91; STAT1
Species	Human
GeneID (Human)	6772
GeneID	6772
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 513-712 of human STAT1 (NP_009330.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42224
Protein Weight	87kDa
Shipping	Ice bag