

Anti-Stat2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P52630

Item Number: CM0017918



Introduction

Rabbit polyclonal antibody targeting human STAT2, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for interferon signaling and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	STAT2 Rabbit pAb
Remarks/Alias	P113; IMD44; ISGF-3; PTORCH3; STAT113; STAT2
Species	Human
GeneID (Human)	6773
GeneID	6773
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-706 of human STAT2 (NP_005410.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52630
Protein Weight	98kDa
Shipping	Ice bag