

## Anti-Stat3 Rabbit Monoclonal Antibody For Wb And Elisa - P40763

Item Number: CM0003000



### Introduction

High-quality rabbit monoclonal antibody against human STAT3, validated for Western blot and ELISA. Detects human and mouse STAT3, molecular weight 88 kDa. Ideal for JAK/STAT signaling studies, inflammation, and cell cycle research.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] STAT3 Rabbit mAb                                                                                       |
| Remarks/Alias    | APRF; HIES; ADMIO; ADMIO1; T3                                                                                         |
| Species          | Human                                                                                                                 |
| Gene ID (Human)  | 6774                                                                                                                  |
| Gene ID          | 6774                                                                                                                  |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 600-700 of human STAT3 (P40763) |
| Source           | Rabbit                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse                                                                                                          |
| SWISS            | P40763                                                                                                                |
| Protein Weight   | 88kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |