

Anti-Stat3 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P40763

Item Number: CM0005270



Introduction

Rabbit monoclonal STAT3 antibody validated for WB, IF/ICC, IHC-P, ELISA. Detects human, mouse, rat, monkey STAT3 (~88kDa). Ideal for JAK/STAT signaling, inflammation, and oncology research.

[Learn More](#)

Parameter	Value
Product Name	STAT3 Rabbit mAb
Remarks/Alias	APRF; HIES; ADMIO; ADMIO1
Species	Human
GeneID (Human)	6774
GeneID	6774
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 611-710 of human STAT3 (NP_003141.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P40763
Protein Weight	88kDa
Shipping	Ice bag