

Anti-Stat3 Rabbit Polyclonal Antibody For Wb And Elisa - P40763

Item Number: CM0020777



Introduction

Rabbit polyclonal antibody targeting human STAT3 (P40763), validated for Western blot and ELISA. Cross-reacts with mouse and rat STAT3. Ideal for studying JAK/STAT signaling, Th17 differentiation, and transcriptional regulation in inflammation and oncology.

[Learn More](#)

Parameter	Value
Product Name	STAT3 Rabbit pAb
Remarks/Alias	APRF; HIES; ADMIO; ADMIO1; STAT3
Species	Human
Gene ID	6774
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-770 of human STAT3 (NP_644805.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40763
Protein Weight	88kDa
Shipping	Ice bag