

Anti-Stat4 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q14765

Item Number: CM0011263

Introduction

Anti-STAT4 Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets the STAT4 transcriptional regulator, critical for Th1 differentiation and interferon signaling. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	STAT4 Rabbit mAb
Remarks/Alias	SLEB11; STAT4
Species	Human
Gene ID	6775
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 649-748 of human STAT4 (NP_003142.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14765
Protein Weight	86kDa
Shipping	Ice bag