

Anti-Stat4 Rabbit Polyclonal Antibody For Wb - Q14765

Item Number: CM0025077

Introduction

STAT4 Rabbit pAb is a polyclonal antibody raised in rabbit against recombinant human STAT4 protein, validated for WB, IHC-P, IF/ICC, ELISA. Detects Human, Mouse, Rat STAT4. Ideal for studying immune response and T-cell differentiation pathways.

[Learn More](#)

Parameter	Value
Product Name	STAT4 Rabbit pAb
Remarks/Alias	SLEB11; STAT4
Species	Human
GeneID (Human)	6775
GeneID	6775
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 569-748 of human STAT4 (NP_003142.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14765
Protein Weight	86kDa
Shipping	Ice bag