


Anti-Stat5A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P42229

Item Number: CM0026228



Introduction

This STAT5A rabbit polyclonal antibody is validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat STAT5A. Ideal for studying signal transduction and transcription regulation in lactation and cancer.

[Learn More](#)

Parameter	Details
Product Name	STAT5A Rabbit pAb
Remarks/Alias	MGF; STAT5; STAT5A
Species	Human
GeneID (Human)	6776
GeneID	6776
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42229
Protein Weight	91kDa
Shipping	Ice bag