

Anti-Stat6 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P42226

Item Number: CM0002234



Introduction

High-specificity anti-STAT6 rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Detects human and mouse STAT6, a key transcription factor in IL-4/IL-3 signaling, supporting cytokine pathway research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	STAT6 Rabbit mAb
Remarks/Alias	STAT6B; STAT6C; D12S1644; IL-4-STAT; STAT6
Species	Human
GeneID (human)	6778
GeneID	6778
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 747-847 of human STAT6 (P42226-1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42226
Protein Weight	94kDa
Shipping	Ice bag