

Anti-Stat6 Rabbit Polyclonal Antibody For Wb, Elisa - P42226

Item Number: CM0019350



Introduction

High-quality rabbit polyclonal antibody against human STAT6, validated for Western blot and ELISA. Ideal for research in IL-4/IL-3 signaling, transcription regulation, and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	STAT6 Rabbit pAb
Remarks/Alias	STAT6; D12S1644; IL-4-STAT; STAT6B; STAT6C; signal transducer and activator of transcription 6
Species	Human
Gene ID (Human)	6778
Gene ID	6778
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 748-847 of human STAT6(NP_003144.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42226
Protein Weight	74kDa/81kDa/94kDa
Shipping	Ice bag