

Anti-Stat6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P42226

Item Number: CM0015383



Introduction

Rabbit polyclonal antibody targeting STAT6, validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat STAT6 for studies on IL-4/STAT6 signaling and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	STAT6 Rabbit pAb
Remarks/Alias	STAT6B; STAT6C; D12S1644; IL-4-STAT; STAT6
Species	Human
GeneID (Human)	6778
GeneID	6778
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42226
Protein Weight	94kDa
Shipping	Ice bag