

Anti-Stathmin 1 Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P16949

Item Number: CM0011585

Introduction

Rabbit monoclonal antibody targeting human Stathmin 1 (OP18, pp17) for WB, IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for microtubule dynamics, leukemia, and neurogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Stathmin 1 Rabbit mAb
Remarks/Alias	Lag; SMN; OP18; PP17; PP19; PR22; LAP18; C1orf215; Stathmin 1
Species	Human
GeneID (Human)	3925
GeneID	3925
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 70-149 of human Stathmin 1 (P16949).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16949
Protein Weight	17kDa
Shipping	Ice bag