

Anti-Stathmin 1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P16949

Item Number: CM0022108

Introduction

This rabbit polyclonal antibody targets Stathmin 1 (Op18/p19) for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Stathmin 1 regulates microtubule dynamics and is overexpressed in acute leukemia. Ideal for cancer research.

[Learn More](#)

Parameter	Specification
Product Name	Stathmin 1 Rabbit pAb
Remarks/Alias	Lag; SMN; OP18; PP17; PP19; PR22; LAP18; C1orf215; Stathmin 1
Species	Human
GeneID (Human)	3925
GeneID	3925
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Stathmin 1 (NP_005554.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16949
Protein Weight	17kDa
Shipping	Ice bag