

Anti-Stathmin 1 Rabbit Polyclonal Antibody For Wb And Elisa - P16949

Item Number: CM0015757

Introduction

High-quality rabbit polyclonal antibody against human Stathmin 1. Validated for Western blot and ELISA, cross-reacts with mouse. Suitable for research on microtubule dynamics and cancer.

[Learn More](#)

Parameter	Value
Product Name	Stathmin 1 Rabbit pAb
Remarks/Alias	Lag, SMN, OP18, PP17, PP19, PR22, LAP18, C1orf215, Stathmin 1
Species	Human
Gene ID	3925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-149 of human Stathmin 1 (NP_981944.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16949
Protein Weight	17kDa
Shipping	Ice bag