

Anti-Stau2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9NuI3

Item Number: CM0029605

Introduction

Rabbit polyclonal antibody targeting human STAU2, a double-stranded RNA-binding protein involved in mRNA transport to dendrites. Validated for WB, IF/ICC, ELISA in human, mouse, and rat samples. Ideal for neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	STAU2 Rabbit pAb
Remarks/Alias	39K2; 39K3; STAU2
Species	Human
GeneID (Human)	27067
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human STAU2 (NP_001157852.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NUL3
Protein Weight	63kDa
Shipping	Ice bag