

Anti-Staufen Polyclonal Antibody For Wb, Ihc-P, Elisa - 095793

Item Number: CM0026519



Introduction

Rabbit polyclonal antibody against human Staufen (STAU1), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Targets double-stranded RNA-binding protein involved in mRNA localization and viral particle production. Ideal for gene expression and host-virus interaction research.

[Learn More](#)

Parameter	Value
Product Name	Staufen Rabbit pAb
Remarks/Alias	STAU; PPP1R150; Staufen
Species	Human
GeneID (Human)	6780
GeneID	6780
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 426-577 of human Staufen (NP_059347.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95793
Protein Weight	63kDa
Shipping	Ice bag