

Anti-Stc1 Rabbit Polyclonal Antibody For Wb And Elisa - P52823

Item Number: CM0022828

Introduction

Rabbit polyclonal antibody targeting human STC1 (Stanniocalcin-1). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for studying renal phosphate regulation and hypercalcemia pathways.

[Learn More](#)

Parameter	Value
Product Name	STC1 Rabbit pAb
Remarks/Alias	STC; STC1
Species	Human
Gene ID (Human)	6781
Gene ID	6781
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P52823
Protein Weight	28kDa
Shipping	Ice bag