

Anti-Stc1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P52823

Item Number: CM0026939

Introduction

Rabbit polyclonal antibody targeting human STC1 (Stanniocalcin-1), validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Useful for studies of renal phosphate reabsorption and hormone signaling.

[Learn More](#)

Parameter	Value
Product Name	STC1 Rabbit pAb
Remarks/Alias	STC; STC1
Species	Human
GeneID (Human)	6781
GeneID	6781
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-247 of human STC1 (NP_003146.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52823
Protein Weight	28kDa
Shipping	Ice bag