

Anti-Steap1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uhe8

Item Number: CM0021399



Introduction

Rabbit polyclonal antibody targeting STEAP1 (six-transmembrane epithelial antigen of prostate 1). Suitable for Western blot and ELISA; reacts with human, mouse, and rat. Useful for prostate tumor and iron metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	STEAP1 Rabbit pAb
Remarks/Alias	STEAP; PRSS24; STEAP1
Species	Human
GeneID (Human)	26872
GeneID	26872
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human STEAP1 (NP_036581.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHE8
Protein Weight	40kDa
Shipping	Ice bag