

Anti-Steroidogenic Factor-1 (Sf-1/Nr5A1) Rabbit Monoclonal Antibody For Wb, Elisa - Q13285

Item Number: CM0003178



Introduction

Steroidogenic factor-1 (SF-1/NR5A1) rabbit monoclonal antibody for Western blot (WB) and ELISA applications. Reacts with human and mouse. Targets SF-1 (52kDa), a transcription factor essential for sexual differentiation and steroidogenic tissue formation.

[Learn More](#)

Parameter	Value
Product Name	Steroidogenic factor-1 (SF-1/NR5A1) Rabbit mAb
Remarks/Alias	ELP; SF1; FTZ1; POF7; SF-1; AD4BP; FTZF1; SPGF8; SRXX4; SRXY3; hSF-1; Steroidogenic factor-1 (SF-1/NR5A1)
Species	Human
GeneID (Human)	2516
GeneID	2516
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 351-461 of human Steroidogenic factor-1 (SF-1/NR5A1) (Q13285)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13285
Protein Weight	52kDa
Shipping	Ice bag