

Anti-Steroidogenic Factor-1 (Sf-1/Nr5A1) Rabbit Polyclonal Antibody For Wb & Elisa - Q13285

Item Number: CM0021711



Introduction

Rabbit polyclonal antibody targeting human Steroidogenic factor-1 (SF-1/NR5A1). Validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Ideal for research on steroidogenic tissue development and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	Steroidogenic factor-1 (SF-1/NR5A1) Rabbit pAb
Remarks/Alias	ELP; SF1; FTZ1; POF7; SF-1; AD4BP; FTZF1; SPGF8; SRXX4; SRXY3; hSF-1; Steroidogenic factor-1 (SF-1/NR5A1)
Species	Human
Gene ID	2516
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Steroidogenic factor-1 (SF-1/Steroidogenic factor-1 (SF-1/NR5A1)) (NP_004950.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13285
Protein Weight	52kDa
Shipping	Ice bag