


Anti-Sterol Carrier Protein 2 (Scp2) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P22307

Item Number: CM0005044



Introduction

Recombinant rabbit monoclonal antibody targeting human Sterol Carrier Protein 2 (SCP2, P22307). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for lipid metabolism and peroxisomal oxidation research.

[Learn More](#)

Parameter	Value
Product Name	Sterol carrier protein 2 Rabbit mAb
Remarks/Alias	NLTP; SCOX; SCPX; SCP-2; SCP-X; NSL-TP; SCP-CHI; Sterol carrier protein 2
Species	Human
Gene ID (Human)	6342
Gene ID	6342
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Sterol carrier protein 2 (P22307)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22307
Protein Weight	59kDa
Shipping	Ice bag