

Anti-Stim1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q13586

Item Number: CM0005439



Introduction

Recombinant anti-STIM1 rabbit monoclonal antibody validated for WB, IHC-P, IP, and ELISA applications. Detects human STIM1, a calcium sensor critical for store-operated calcium entry. Ideal for calcium signaling and ion transport research.

[Learn More](#)

Parameter	Value
Product Name	STIM1 Rabbit mAb
Remarks/Alias	GOK; TAM; TAM1; IMD10; STRMK; D11S4896E; STIM1
Species	Human
GeneID (Human)	6786
GeneID	6786
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 530-685 of human STIM1 (NP_003147.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q13586
Protein Weight	77kDa
Shipping	Ice bag