

Anti-Stim1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q13586

Item Number: CM0031049



Introduction

Rabbit polyclonal antibody targeting human STIM1 (SWISS Q13586), validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for Ca²⁺ signaling and store-operated calcium entry (SOCE) research.

[Learn More](#)

Parameter	Value
Product Name	STIM1 Rabbit pAb
Remarks/Alias	GOK; TAM; TAM1; IMD10; STRMK; D11S4896E; STIM1
Species	Human
GeneID (Human)	6786
GeneID	6786
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 530-685 of human STIM1 (NP_003147.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13586
Protein Weight	77kDa
Shipping	Ice bag