

Anti-Stim1 Rabbit Polymab® For Wb, Ihc-P, Ip, Elisa - Q13586

Item Number: CM0006300



Introduction

Recombinant rabbit monoclonal antibody targeting human STIM1 (aa 530-685) for Western Blot, IHC-P, IP, and ELISA. Designed for researchers studying store-operated calcium entry. Cross-reactive with mouse and rat.

[Learn More](#)

Parameter	Specification
Product Name	STIM1 Rabbit PolymAb®
Remarks/Alias	GOK; TAM; TAM1; IMD10; STRMK; D11S4896E
Species	Human
GeneID (Human)	6786
GeneID	6786
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 530-685 of human STIM1 (NP_003147.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13586
Protein Weight	77kDa
Shipping	Ice bag