

Anti-Stim2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9P246

Item Number: CM0019875



Introduction

Rabbit polyclonal antibody against human STIM2, validated for Western blot, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat STIM2; immunogen corresponds to amino acids 60-180. Suitable for studying store-operated calcium entry.

[Learn More](#)

Parameter	Value
Product Name	STIM2 Rabbit pAb
Remarks/Alias	STIM2
Species	Human
Gene ID	57620
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-180 of human STIM2 (NP_001162588)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P246
Protein Weight	84kDa
Shipping	Ice bag