

Anti-Sting/Tmem173 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q86Wv6

Item Number: CM0005764



Introduction

Rabbit monoclonal antibody against human STING/TMEM173, validated for WB, IHC-P, ELISA. Detects human, mouse, rat. Ideal for innate immunity and autophagy research. Molecular weight 42 kDa.

[Learn More](#)

Parameter	Value
Product Name	STING/TMEM173 Rabbit mAb
Remarks/Alias	ERIS; MITA; MPYS; SAVI; NET23; STING; hMITA; hSTING; TMEM173; STING-beta; STING/TMEM173
Species	Human
GeneID (Human)	340061
GeneID	340061
Immunogen	Recombinant protein of human STING/TMEM173
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86WV6
Protein Weight	42kDa
Shipping	Ice bag