

Anti-Sting/Tmem173 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - Q86Wv6

Item Number: CM0030760



Introduction

Rabbit polyclonal antibody against human STING/TMEM173, validated for WB, IHC-P, IF/ICC, ELISA. Detects STING in human, mouse, and rat samples. Suitable for innate immunity and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	STING/TMEM173 Rabbit pAb
Remarks/Alias	ERIS; MITA; MPYS; SAVI; NET23; STING; hMITA; hSTING; TMEM173; STING-beta; STING/TMEM173
Species	Human
GeneID (Human)	340061
GeneID	340061
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 254-379 of human STING/TMEM173 (NP_938023.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86WV6
Protein Weight	42kDa
Shipping	Ice bag