

Anti-Sting/Tmem173 Rabbit Polyclonal Antibody - Q86Wv6

Item Number: CM0027677



Introduction

High-quality anti-STING/TMEM173 rabbit polyclonal antibody validated for Western blot, immunoprecipitation, immunofluorescence, immunohistochemistry, and ELISA. Cross-reacts with human, mouse, and rat STING, enabling innate immunity and cGAS-STING pathway research.

[Learn More](#)

Parameter	Value
Product Name	STING/TMEM173 Rabbit pAb
Remarks/Alias	ERIS; MITA; MPYS; SAVI; NET23; STING; hMITA; hSTING; TMEM173; STING-beta; STING/TMEM173
Species	Human
GeneID (human)	340061
GeneID	340061
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 139-379 of human STING/TMEM173 (NP_938023.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86WV6
Protein Weight	42 kDa
Shipping	Ice bag