

Anti-Stip1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P31948

Item Number: CM0012993



Introduction

STIP1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, IP, and ELISA. Reacts with human, mouse, rat, monkey. Suitable for stress-induced phosphoprotein 1 co-chaperone research.

[Learn More](#)

Parameter	Value
Product Name	STIP1 Rabbit pAb
Remarks/Alias	HOP; P60; STI1; STI1L; HEL-S-94n; IEF-SSP-3521; STIP1
Species	Human
GeneID (Human)	10963
GeneID	10963
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P31948
Protein Weight	63kDa
Shipping	Ice bag