

Anti-Stip1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P31948

Item Number: CM0003198



Introduction

Rabbit monoclonal antibody against human STIP1 (Hop) validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat STIP1; immunogen: amino acids 2-98 of human STIP1.

[Learn More](#)

Parameter	Value
Product Name	STIP1 Rabbit mAb
Remarks/Alias	HOP; P60; STI1; STI1L; HEL-S-94n; IEF-SSP-3521; STIP1
Species	Human
Gene ID (Human)	10963
Gene ID	10963
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-98 of human STIP1 (P31948)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31948
Protein Weight	63kDa
Shipping	Ice bag