


Anti-Stk19 Polyclonal Antibody For Wb, Ihc-P, Elisa - P49842

Item Number: CM0016473



Introduction

High-quality rabbit polyclonal antibody against STK19/WHR1, validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat STK19, a key factor in transcription-coupled DNA repair.

[Learn More](#)

Parameter	Value
Product Name	STK19 Rabbit pAb
Remarks/Alias	G11; RP1; D6S60; D6S60E; HLA-RP1; STK19
Species	Human
GeneID (Human)	8859
GeneID	8859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 75-364 of human STK19 (NP_004188.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49842
Protein Weight	41kDa
Shipping	Ice bag