


Anti-Stk19 Rabbit Polyclonal Antibody For Wb, Elisa - P49842

Item Number: CM0027611



Introduction

STK19 Rabbit Polyclonal Antibody validated for WB, ELISA in Human, Mouse, Rat. Detects 41kDa STK19/WHR1 winged helix repair factor. Ideal for transcription-coupled repair research.

[Learn More](#)

Parameter	Value
Product Name	STK19 Rabbit pAb
Remarks/Alias	G11; RP1; D6S60; D6S60E; HLA-RP1; STK19
Species	Human
GeneID	8859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 75-364 of human STK19 (NP_004188.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49842
Protein Weight	41kDa
Shipping	Ice bag