

Anti-Stk24 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6E0

Item Number: CM0016873



Introduction

Rabbit polyclonal antibody targeting human STK24 (MST3) for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for apoptosis and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	STK24 Rabbit pAb
Remarks/Alias	MST3; STK3; MST3B; STE20; HEL-S-95; STK24
Species	Human
GeneID (Human)	8428
GeneID	8428
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 312-431 of human STK24 (NP_001027467.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6E0
Protein Weight	49kDa
Shipping	Ice bag