

Anti-Stk26 (Mst4) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9P289

Item Number: CM0023197



Introduction

Rabbit polyclonal antibody against human STK26 (MST4/MASK), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for apoptosis, autophagy, and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	STK26 Rabbit pAb
Remarks/Alias	MASK; MST4; STK26
Species	Human
GeneID (Human)	51765
GeneID	51765
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-416 of human STK26 (NP_057626.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P289
Protein Weight	47kDa
Shipping	Ice bag