

Anti-Stk31 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bxu1

Item Number: CM0018680

Introduction

Anti-STK31 Rabbit Polyclonal Antibody validated for WB and ELISA applications. Recognizes human STK31 (Serine/threonine-protein kinase 31) with cross-reactivity to mouse and rat. Ideal for testis-specific kinase studies.

[Learn More](#)

Parameter	Value
Product Name	STK31 Rabbit pAb
Remarks/Alias	TDRD8; SGK396; STK31
Species	Human
GeneID (Human)	56164
GeneID	56164
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 890-1019 of human STK31 (NP_113602.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BXU1
Protein Weight	116kDa
Shipping	Ice bag