

Anti-Stk33 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Byt3

Item Number: CM0002089



Introduction

Rabbit monoclonal STK33 antibody validated for WB, IF/ICC, ELISA. Detects human, mouse STK33 (~58 kDa), a serine/threonine kinase required for spermatid differentiation and male fertility.

[Learn More](#)

Parameter	Value
Product Name	STK33 Rabbit mAb
Remarks/Alias	STK33
Species	Human
GeneID (Human)	65975
GeneID	65975
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 399-514 of human STK33 (Q9BYT3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BYT3
Protein Weight	58kDa
Shipping	Ice bag