

Anti-Stk36 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nrp7

Item Number: CM0017013

Introduction

Rabbit polyclonal antibody targeting human STK36 (Fused homolog), suitable for WB and ELISA, with cross-reactivity to mouse. Key tool for smoothened pathway, G1I transcription factor regulation, and cilia research.

[Learn More](#)

Parameter	Value
Product Name	STK36 Rabbit pAb
Remarks/Alias	FU, CILD46, STK36
Species	Human
GeneID (Human)	27148
GeneID	27148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-450 of human STK36 (NP_056505.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NRP7
Protein Weight	144kDa
Shipping	Ice bag