


Anti-Stk38 Polyclonal Antibody For Wb, If/Icc, Elisa - Q15208

Item Number: CM0029043



Introduction

Rabbit polyclonal antibody against STK38 (NDR1), tested for WB, IF/ICC, ELISA in human, mouse, rat. Targets a kinase involved in MAPK signaling and DNA damage response. Suitable for cell signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	STK38 Rabbit pAb
Remarks/Alias	NDR; NDR1; STK38
Species	Human
Gene ID	11329
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 366-465 of human STK38 (NP_009202.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15208
Protein Weight	54kDa
Shipping	Ice bag