


Anti-Stk38L Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2H1

Item Number: CM0017534



Introduction

High-quality rabbit polyclonal antibody targeting STK38L (NDR2) for WB, IF/ICC, and ELISA.

Validated in human, mouse, and rat samples. Ideal for studying neuronal differentiation and kinase signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	STK38L Rabbit pAb
Remarks/Alias	NDR2; STK38L
Species	Human
Gene ID (Human)	23012
Gene ID	23012
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 165-464 of human STK38L (NP_055815.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2H1
Protein Weight	54kDa
Shipping	Ice bag